

*Release the Inner Goddess*

QUESTIONS OR COMMENTS  
For Replacement Parts or Facial Scrub



ARTEMIS  
WOMAN™

ARTEMIS WOMAN, LLC  
P. O. Box 1037, Weston CT 06883  
1 866 757 7776  
[www.artemiswoman.com](http://www.artemiswoman.com)

Unit MADE IN CHINA  
Facial Scrub MADE IN USA

## MICRODERMABRASION SYSTEM WITH PORE CLEANSING



 INSTRUCTION BOOKLET

FOR A REVOLUTIONARY NEW LINE OF  
SKIN CARE PRODUCTS BY



ARTEMIS  
WOMAN™

HC-100 ♦ HG-1000

### HOW TO INSERT BATTERY

Press both sides of the battery cover located at the lower end of the body and pull out. Insert two “AA” batteries in the direction indicated in the battery compartment and replace the cover.

### TO REMOVE AND INSTALL THE ATTACHMENTS

Select the attachment you desire and snap tightly in the unit. To remove simply pull out the attachment.

### TO SELECT THE PROPER ACTION

The action selector enables you to select either Exfoliate or Pore Cleanse when you use the attachments. Turn action selector to Exfoliate for use with Cleansing Brush, Massaging Roller and Sponge Applicator. While you use the Pore Cleansing Cups, turn action selector to Pore Cleanse.



### CLEANSING BRUSH

First turn action selector to Rotation before turning on the unit. Place Cleansing Brush on microdermabrasion unit. The Cleansing Brush will help remove facial dirt, dead surface cells and simultaneously activate facial cells.

Apply facial scrub to the Cleansing Brush and place brush gently on your face. Turn unit on, select desired speed, and gently move brush with a light circular motion. Avoid contact with the delicate areas around the eyes. Turn unit off, remove brush from unit, and rinse the brush with clean water.



Exfoliating  
Brush



Sponge Applicator

### SPONGE APPLICATOR

If you prefer, you may cleanse with the Sponge applicator. Place Sponge Applicator on unit. Apply facial scrub, turn unit on, select desired speed, and the soft applicator will help you to achieve a clean, smooth, even and perfect finish to your skin.



## PORE CLEANSER

Remember to turn action selector to Pore Cleanse before you use the Pore Cleanser. There are two kinds of attachments provided. The larger pore cleanser is intended for the face while the smaller one is intended for the nose. The pore cleanser will help loosen surface impurities. To avoid bruising, do not use the Pore Cleanse switch without the Pore Cleansing attachment.



Small Pore  
Cleanser



Large Pore  
Cleanser



Facial Massager

## MASSAGING ROLLER

Apply your regular toner and moisturizer over the face and neck. Place Massaging Roller on unit. Turn unit on, select desired speed, and gently apply the Massaging Roller over face to relieve facial tension and to maximize the effects of the moisturizing treatment. Turn unit off.



## CLEANING, RINSING, AND DRYING

1. Do not submerge appliance in water. Remove attachments prior to cleaning.
2. Always rinse the attachments thoroughly in clean water after use.
3. Detach the Pore Cleanser attachment to clean the inside cavity with clean water and replace the cover.
4. Make sure to dry every attachment with a dry clean cloth after cleansing.
5. Use damp cloth to clean the main power unit.
6. Do not use benzene, thinners, or other organic solvents to clean, because of the risk of damaging the unit.



## CAUTIONS

1. Do not stay on any one area of the face for any extended period of time with any one of the attachments.
2. Do not press the skin too hard as it will slow down the rotation of the unit and could possibly damage your skin.
3. Do not drop or strike the main power unit against any hard object as it will create scratches and damage the unit.
4. Never wash or immerse the main power unit in water or any other liquids, for it is not waterproof and could result in a short circuit.
5. Do not leave a dead battery in your main power unit for a long time as it may leak and corrode the internal contacts.
6. Remember to switch off the unit before you change the attachments.



## WARNINGS

1. Anyone suffering from acne or any other skin condition should consult a doctor before using the facial microdermabrasion product.
2. Avoid contact with the delicate areas around the eyes.
3. Do not hold the Pore Cleanser in one area for an extended period of time. We suggest no more than one minute on the nose and no more than three to four minutes on the cheeks, chin, and forehead. Discoloration or bruising could appear temporarily on your skin, especially if held against the skin for a prolonged period of time. Use the Small Pore Cleansing Cup only on the nose area.
4. The facial scrub is intended for external use only. Do not ingest.

